**KLINION®**

**L-Mesitran®**

**Vernieuwd**

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### Referenties:

- 1.: J Wound Care/Why do some cavity wounds treated with honey heal without scarring. Vol 11. (2) pg 52. 2002
- 2.: J Wound Care/The role of honey in the management of wounds. Vol 8. (8) pg 415. 1999
- 3.: J Roy Soc Med/Antibacterial activity of honey against Staphylococcus from infected wounds. Vol 92. pg 283. 1999
- 4.: J Postgard Med/Effect of honey on multidrug resistant organisms and its synergistic action with three common antibiotics. Vol 44. pg 93. 1998
- 5.: J Obstet Gyn/Radical operation for carcinoma of the vulva. Vol 77. (11) pg 1037. 1970
- 6.: J Surgery/Recent advances in the management of Fournier's gangrene. Vol 113. (2) pg 200. 1993
- 7.: J Urology/Genital Fournier's gangrene. Vol 47. (5) pg 734. 1996
- 8.: Brit. J Surgery/Clinical observations on the wound healing properties of honey. Vol 75. pg 679. 1988
- 9.: Brit. J Surgery/Topical application of honey in treatment of burns. Vol 78. (4) pg 497. 1991
- 10.: Brit. J Plast. Surg/Honey impr. gauze versus polyurethane film in the treatment of burns. Vol 46. (4) pg 322. 1993
- 11.: J Orthop Surgery/Clinical and bacteriological results in wounds treated with honey. Vol 7. (2) pg 202. 1993
- 12.: J Burns/A clinical and histological study of superficial burn wound healing with honey. Vol 24. (2) pg 157. 1998
- 13.: J Obstet Gyn/Speeding up the healing of burns with honey. Vol 32. pg 381. 1992
- 14.: J Trop Med Hyg/The pattern and outcome of burn injuries, a review of 156 cases. Vol 97. pg 108. 1994
- 15.: J Burn Care/The efficacy of honey in inhibiting strains of Pseudomonas aer. from infected burns. Vol 23. (6) pg 366-370. 2002
- 16.: J Appl Microbiol/The sensitivity to honey of Gram-positive cocci of clinical significance isolated from wounds. Vol 93. (5) pg 857-863. 2002
- 17.: J Dermatology Treat/Healing of an MRSA-colonized, hydroxyurea-induced leg ulcer with honey. Vol 12. (1) pag 33-36. 2001
- 18.: Acta Paediatr/Local application of honey for treatment of neonatal postoperative wound infection. Vol 87. (4) pg 429-432. 1998
- 19.: J R Soc Med/Antibacterial activity of honey against strains of Staphylococcus aureus from infected wounds. Vol 92. (6) pg 283-285. 1999
- 20.: Austr. Wound Man. Ass. Conf./Honey has potential as a dressing for wounds infected with MRSA, March 1998
- 21.: J Wound Care/The use of honey as an antiseptic in managing Pseudomonas infection. Vol 8. (4) pg 161-164. 1999
- 22.: J Pharm Pharmacol/The effect of gamma-irradiation on the antibacterial activity of honey. Vol 48. pg 1206. 1996
- 23.: Ostomy Wound Manage/Re-introducing Honey in the Management of Wounds and Ulcers. Vol 48. (11) pg 28-40. 2002
- 24.: Am. J Clin Dermatol/Potential of honey in the treatment of wounds and burns. Vol 2. (1) pg 13-19. 2001
- 25.: Ligament/Honing en wondzorg back to the future. Vol 33. (2) pg 17. 2002
- 26.: Austr. J. of Wound Man./The evidence for honey promoting wound healing Vol 6. (4) pg 148-158. 1998
- 27.: Nurs Stand/The use of honey in wound management. Vol 15. (11) pg 63. 2000
- 28.: Arch Surgery/Protective covering of surgical wounds with honey impedes Tumor Implantation. vol 135. 2000
- 29.: Supplement/Leids Cytologisch en Pathologisch Laboratorium/Boon ME/Onderzoek naar effect van honingzalf voor vaginaal gebruik. Vol 26. febr. 2004
- 30.: Nagelschimmels in de Leidse praktijk/ Leids Cytologisch en Pathologisch Laboratorium/Boon M.E. Middag-Broekman J.H.FF. Jaarverslag LCPL 2003
- 31.: Addition of antioxidants and Polyethylene Glycol 4000 enhances the healing property of honey in burns. Subrahmanyam N. Department of Surgery, Dr. Vaishampayan Memorial Medical College, India Annals of Burns and Fire Disastus, vol 9, nr 2, juni 1996.

**Studies over honing/L-Mesitran m.b.t. de effectiviteit in wondgenezing en de antibacteriële werking. (abstracts National Library of Medicine PubMed)**

**Topical application of honey in treatment of burns.**

Br J Surg 1991; 78(4):497-498. Subrahmanyam M.

Honey was compared with silver sulfadiazine impregnated gauze (Flamazine) for efficacy as a dressing for superficial burn injury in a prospective randomised controlled trial that was carried out with a total of 104 patients. In the 52 patients treated with honey, 91% of the wounds were rendered sterile within 7 days. In the 52 patients treated with silver sulfadiazine, 7% showed control of infection within 7 days. Healthy granulation tissue was observed earlier in patients treated with honey (means 7.4 versus 13.4 days). The time taken for healing was significantly shorter with the honeytreated group (p<0.001): of the wounds treated with honey 87% healed within 15 days compared with 10% of those treated with silver sulfadiazine. Better relief of pain, less exudation, less irritation of the wound, and a lower incidence of hypertrophic scar and post-burn contracture were noted with the honey treatment. The honey treatment also gave acceleration of epithelialisation at 6 - 9 days, a chemical debridement effect and removal of offensive smell.

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**Healing of an MRSA-colonized, hydroxyurea-induced leg ulcer with honey.**

J Dermatology Treat 2001 March;12(1): 33-6/Natarajan S, Williamson D, Grey J, Harding KG, Cooper RA. Wound Healing Research Unit, University of Wales College of Medicine, Heath Park, Cardiff, UK.

BACKGROUND: With the everincreasing emergence of antibiotic-resistant pathogens, in particular methicillin-resistant Staphylococcus aureus (MRSA) in leg ulcers, a means of reducing the bacterial bioburden of such ulcers, other than by the use of either topical or systemic antibiotics, is urgently required. METHODS: We report the case of an immunosuppressed patient who developed a hydroxyurea-induced leg ulcer with subclinical MRSA infection which was subsequently treated with topical application of manuka honey, without cessation of hydroxyurea or cyclosporin. RESULTS: MRSA was eradicated from the ulcer and rapid healing was successfully achieved. CONCLUSION: Honey is recognized to have anti-bacterial properties, and can also promote effective wound healing. A traditional therapy, therefore, appears to have enormous potential in solving new problems.

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**Mesitran clinical reports study 2002**

Drs. Jan Vandeputte Clinical Nursing Consulting International bvba

Three separate studies (139 patients) were used to obtain some relevant data concerning the effect of Mesitran, (a honey based ointment) in acute and chronic wounds. The studies were executed in different facilities, 2 general hospitals, 12 nursing homes and home care facilities. The wounds that were evaluated were acute wounds (bruises and burns) and chronic wounds (diabetic, venous insufficiencies, arterial diseases and pressure ulcers). The best results were seen in small 2n degree burns (47% faster) and venous ulcers (42% faster). Mesitran was able to heal skin tears 29% quicker than any other product. Less difference (14%) was seen in the healing of pressure sores. The Mesitran, wound care ointment has a great potential to heal acute and chronic wounds.

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**Honey has potential as a dressing for wounds infected with MRSA,**

The Second Australian Wound Management Association Conference, Brisbane, Australia, 18-21 March, 1998, 1998. Molan P, Brett M.

Honey provides a moist healing environment yet prevents bacterial growth even when wounds are heavily infected. It is a very effective means of quickly rendering heavily infected wounds sterile, without the side-effects of antibiotics, and it is effective against antibiotic-resistant strains of bacteria.

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**Re-introducing Honey in the Management of Wounds and Ulcers - Theory and Practice.**

Ostomy Wound Manage 2002 Nov; 48(11):28-40/Molan PC. Department of Biological Sciences, University of Waikato, Private Bag 3105, Hamilton, New Zealand.

Dressing wounds with honey, a standard practice in past times, went out of fashion when antibiotics came into use. Because antibiotic-resistant bacteria have become a widespread clinical problem, a renaissance in honey use has occurred. Laboratory studies and clinical trials have shown that honey is an effective broad-spectrum antibacterial agent that has no adverse effects on wound tissues. As well as having an antibacterial action, honey also provides rapid autolytic debridement, deodorizes wounds, and stimulates the growth of wound tissues to hasten healing and start the healing process in dormant wounds. Its anti-inflammatory activity rapidly reduces pain, edema, and exudate and minimizes hypertrophic scarring. It also provides a moist healing environment for wound tissues with no risk of maceration of surrounding skin and completely prevents adherence of dressings to the wound bed so no pain or tissue damage is associated with dressing changes. Using appropriate dressing

practice overcomes potential messiness and handling problems.

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**The efficacy of honey in inhibiting strains of Pseudomonas aeruginosa from infected burns.**

J Burn Care Rehabil 2002 Nov-Dec; 23(6):366-70/ Cooper RA, Halas E, Molan PC. Centre for Biomedical Sciences, School of Applied Sciences, University of Wales Institute, Cardiff, Wales, United Kingdom.

Because there is no ideal therapy for burns infected with Pseudomonas aeruginosa, there is sufficient need to investigate the efficacy of alternative anti-pseudomonal interventions. Honey is an ancient wound remedy for which there is modern evidence of efficacy in the treatment of burn wounds, but limited evidence for the effectiveness of its antibacterial activity against Pseudomonas. We tested the sensitivity of 17 strains of P. aeruginosa isolated from infected burns to two honeys with different types of antibacterial activity, a pasture honey and a manuka honey, both with median levels of activity. All strains showed similar sensitivity to honey with minimum inhibitory concentrations below 10% (vol/vol); both

honeys maintained bactericidal activity when diluted more than 10-fold. Honey with proven antibacterial activity has the potential to be an effective treatment option for burns infected or at risk of infection with P. aeruginosa.

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**The sensitivity to honey of Gram-positive cocci of clinical significance isolated from wounds.**

J Appl Microbiol 2002 Nov; 93(5):857- 863/ Cooper RA, Molan PC, Harding KG. Centre for Biomedical Sciences, School of Applied Sciences, University of Wales Institute Cardiff, Llandaff Campus, Cardiff, Wales.

AIMS: To determine the sensitivity to honey of Gram-positive cocci of clinical significance in wounds and demonstrate that inhibition is not exclusively due to osmotic effects. METHODS AND RESULTS: Eighteen strains of methicillin-resistant Staphylococcus aureus (MRSA) and seven strains of vancomycin-sensitive enterococci (VRE) were isolated from infected wounds and 20 strains of vancomycin-resistant enterococci were isolated from hospital environmental surfaces. Using an agar incorporation technique to determine the minimum inhibitory concentration (MIC), their sensitivity to two natural honeys of median levels of antibacterial activity was established and compared with an artificial honey solution. For all of the strains tested, the MIC values against manuka and pasture honey were below 10% (v/v), but concentrations of artificial honey at least three times higher were required to achieve equivalent inhibition in vitro. Comparison of the MIC values of antibiotic-sensitive strains with

their respective antibiotic-resistant strains demonstrated no marked differences in their susceptibilities to honey.

CONCLUSIONS: The inhibition of bacteria by honey is not exclusively due to osmolality. For the Gram-positive cocci tested, antibiotic-sensitive and -resistant strains showed similar sensitivity to honey. SIGNIFICANCE AND IMPACT OF THE STUDY: A possible role for honey in the treatment of wounds colonized by antibiotic-resistant bacteria is indicated.

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**Local application of honey for treatment of neonatal postoperative wound infection.**

Acta Paediatr 1998 April; 87(4): 429- 32/Vardi A, Barzilay Z, Linder N, Cohen HA, Paret G, Barzilay A. Paediatric Intensive Care Unit, Chaim Sheba Medical Centre, Tel Hashomer, Israel.

Honey has been described in ancient and modern medicine as being effective in the healing of various infected wounds. In this report we present our experience in nine infants with large, open, infected wounds that failed to heal with conventional treatment. Conventional treatment was defined as having failed if after > or = 14 d of intravenous antibiotic and cleaning the wound with chlorhexidine 0.05% W/V in aqueous solution and fusidic acid ointment the wound was still open, oozing pus, and swab cultures were positive. All infants showed marked clinical improvement after 5 d of treatment with topical application of 5-10 ml of fresh unprocessed honey twice daily. The wounds were closed, clean and sterile in all infants after 21 d of honey appli-

cation. There were no adverse reactions to the treatment.

We conclude that honey is useful in the treatment of post-surgical wounds that are infected and do not respond to conventional systemic and local antibiotic treatment.

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**Antibacterial activity of honey against strains of Staphylococcus aureus from infected wounds.**

J R Soc Med 1999 June; 92(6): 283-5/ Cooper RA, Molan PC, Harding KG. School of Biomedical Sciences, University of Wales Institute, Cardiff, UK.

The antibacterial action of honey in infected wounds does not depend wholly on its high osmolality. We tested the sensitivity of 58 strains of coagulase-positive Staphylococcus aureus, isolated from infected wounds, to a pasture honey and a manuka honey. There was little variation between the isolates in their sensitivity to honey: minimum inhibitory concentrations were all between 2 and 3% (v/v) for the manuka honey and between 3 and 4% for the pasture honey. Thus, these honeys would prevent growth of S. aureus if diluted by body fluids a further seven-fold to fourteen-fold beyond the point where their osmolality ceased to be completely inhibitory. The antibacterial action of the pasture honey relied on release of hydrogen peroxide, which in vivo might be reduced by catalase activity in tissues or blood. The action of manuka honey stems partly from a phytochemical component, so this type of honey might be more effective in vivo. Comparative clinical trials with standardized honeys are needed.

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**The effect of gamma-irradiation on the antibacterial activity of honey.**

J Pharm Pharmacol 1996; 48:1206-1209. Molan PC, Allen KL.

**The sterilization of honey with cobalt 60 gamma radiation: a study of honey spiked with Clostridium botulinum and Bacillus subtilis.** Experientia (Basel) 1995; 51:986-989. Postmes T, van den Bogaard AE, Hazen M.

Honey sometimes contains spores of clostridia, which poses a small risk of wound botulism. This objection can be overcome by the use of honey that has been treated by gamma-irradiation, which kills clostridial spores in honey without loss of any of the antibacterial activity. (Mesitran honey is gamma-irradiated)

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**Potential of honey in the treatment of wounds and burns.**

Am J Clin Dermatol 2001;2(1): 13-9/Molan PC. Honey Research Unit, University of Waikato, Hamilton, New Zealand. pmolan@waikato.ac.nz

There has been a renaissance in recent times in the use of honey, an ancient and traditional wound dressing, for the treatment of wounds, burns, and skin ulcers. In the past decade there have been many reports of case

studies, experiments using animal models, and randomized controlled clinical trials that provide a large body of very convincing evidence for its effectiveness, and biomedical research that explains how honey produces such good results. As a dressing on wounds, honey provides a moist healing environment, rapidly clears infection, deodorizes, and reduces inflammation, edema, and exudation. Also, it increases the rate of healing by stimulation of angiogenesis, granulation, and epithelialization, making skin grafting unnecessary and giving excellent cosmetic results.

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**The evidence for honey promoting wound healing**  
**A brief review of the use of honey as a clinical dressing.**

Molan, P. C. Primary Intention (The Australian Journal of Wound Management) 6 (4) 148-158 (1998).

The clinical observations recorded are that infection is rapidly cleared, inflammation, swelling and pain are quickly reduced, odour is reduced, sloughing of necrotic tissue is induced, granulation and epithelialisation are hastened, and healing occurs rapidly with minimal scarring.

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**Inflammation.**

Ryan GB, Majno G. Kalamazoo, Michigan: Upjohn, 1977:80.

It provides a supply of glucose for leucocytes, essential for the 'respiratory burst' that produces hydrogen peroxide, the dominant component of the antibacterial activity of macrophages. The acidity of honey (typically below pH 4) may also assist in the antibacterial action of macrophages, as an acid pH inside the vacuole is involved in killing ingested bacteria.

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**The use of honey in wound management.**

Nurs Stand 2000 Nov 29- Dec 5; 15(11):63-8/Dunford C, Cooper R, Molan P, White R. Salisbury District Hospital.

Honey has been used as a wound treatment for more than 2,000 years. Greater scientific understanding of how it works, particularly as an antibacterial agent, has led practitioners to reconsider the therapeutic value of honey. Once honey is commercially available as a regulated product in the UK, practitioners will have access to an effective, alternative wound treatment. Specific, sterilised honeys intended for wound care will provide a safe natural product to manage colonised or infected wounds that would otherwise remain unresponsive to treatment.

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**Why do some cavity wounds treated with honey or sugar paste heal without scarring?**

J. Topham, MPhil, BPharm, MRPharmS, Journal of Wound Care, vol 11, no 2, February 2002.

Saccharides at the wound surface encourage the production of hyaluronic acid from glucose, simultaneously suppressing the formation of fibroforming collagens. The sugar preparations at the wound bed create an environment that enables woundhealing proteoglycans to exert their effects without producing excessive quantities of collagens. The saccharide attachments to the nascent collagen may result in branching of the triplestranded helical structure of collagens. This will pro-

duce the meshlike scaffold structure of collagen type IV. Vitamin C is important for the prevention of scars when sugar or honey are applied to wounds.

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**Protective covering of surgical wounds with honey impedes Tumor Implantation.**

Arch Surgery, vol 135, december 2000. Ismail Hamzaoglu et al. Departments of Surgery, Pathology and Medical Biology, University of Istanbul, Turkey.

Tumor implantation (TI) development at the surgical wound following cancer surgery is still an unresolved concern. Trocar site recurrence (TSR), which is likely a form of TI, has become one of the most controversial topics and, with the widespread acceptance of laparoscopic surgery, has caused renewed interest in questions about IT. Physiological and chemical properties of honey might prevent Tumor Implantation when applied locally.

Honey could be used on trocar wounds of patients with malignant disease to prevent Trocar Site Recurrence (TSR). It also may prevent infection with its bactericidal effects. Wound healing also may even be enhanced with topical application of honey. Honey may also provide benefit in conventional oncological surgery where TI is predictable.

Preventive covering of surgical wounds with honey seems to be a harmless procedure and may constitute at least a partial barrier that might overcome Tumor cell invasion.

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**The use of honey as an antiseptic in managing Pseudomonas infection.**

J Wound Care 1999 April; 8(4):161-4/Cooper R, Molan P. School of Biomedical Sciences, University of Wales Institute, Cardiff, UK.

A laboratory study was undertaken to extend existing knowledge about the effectiveness of the antibacterial properties of honey against pseudomonas. To date, sensitivity testing has used non-standardised honeys, which may vary greatly in their antibacterial potency. Pure cultures of Pseudomonas spp, isolated from swabs from 20 infected wounds, were inoculated on the surface of nutrient agar plates containing various concentrations of honey in the medium. Two types of honey were used, a manuka honey and a pasture honey, each selected to have antibacterial activity close to the median for each type. The minimum inhibitory concentration of the manuka honey for the 20 isolates ranged from 5.5-8.7% (v/v) (mean 6.9% (v/v), standard deviation 1.3). The minimum inhibitory concentration of the pasture honey for the 20 isolates ranged from 5.8-9.0% (v/v) (mean 7.1% (v/v), standard deviation 1.0). Honeys with an average level of antibacterial activity could be expected to be effective in preventing the growth of pseudomonas on the surface of a wound even if the honey were diluted more than tenfold by exudation from the wound.

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**Addition of antioxidants and polyethylene glycol 4000 enhances the healing property of honey in burns**

Subrahmanyam N. Department of Surgery, Dr. Vaishampayan Memorial Medical College, India

A prospective study was carried out in two groups of randomly allocated 84 patients with partial-thickness burns to compare the effect of honey alone with that of honey fortified with antioxidants (vitamin C and vitamin E) and polyethylene glycol. Burns treated with honey plus antioxidants and PEG, healed earlier than those treated with honey alone. In this study the addition of antioxidants and PEG accelerated the healing properties of honey in burns. In patients treated with honey alone, healthy granulation appeared in the majority of cases after about one week: wounds healed by day 10 days in 20 patients and by day 15 in another 12 patients. Thus 32 patients (76. 1 %) achieved wound healing by day 15. In patients treated with honey PEG/C/E, granulation appeared by day 4 in the majority of cases, with wounds healing by day 6 in 34 patients and by day 9 in another 8 patients. Thus 42 patients (100%) achieved wound healing within 9 days. The time taken for wound healing differed significantly between the two groups. Allergy or other side-effects were not observed in any patient in either group.

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Medeco B.V., Postbus 1555, 3260 BB Oud-Beijerland